

Biax Experiment (rev. 27 June 2019)

Exp. Name: P5356 W6 Frac 5-20

Date/Time: 10/2/19

Operator: Wood/Manogharan

Hydraulics start: _____

Example name: PXXXXBttMatNN

Hydraulics end: _____

Sample Block Thickness w/ no gouge:

___ Steel 5x5 cm, _____ mm

___ Vessel (Small Single Direct)-Frits: _____

___ Titanium 5x5 cm, _____ mm

___ Vessel (Large Single Direct)

___ Steel 10x10 cm, _____ mm

___ Vessel (5x5 Grooved)-Frits: _____

___ Titanium 10x10 cm, _____ mm

Vessel Side Blocks: _____ Empty Block + frits: _____

For Current Calibrations see: ~gps/group/cjm38/default/Calibrations/


Layer Thickness (total on bench): _____ mm Under Load: _____ mm@sample

Material (Qtz, Granite, ?): W6, Pre-fracture (from P5346 intact)

Particle Size, Size Distribution: _____


Load cells:

Contact area: 0.002233036 m²

Load cell name	Calibrations (mV/kN)		Target stress (MPa)	Init. Voltage	Volt. @ load
62 mm H	LG: 18.561	HG: 172.1	Hor: <u>5, 10, 15, 20</u>	10.0000 <u>-0.038339</u>	
<u>44 mm H</u>	LG: 12.3	<u>HG: 123.9</u>			
22 mm H	Gain: 773.6		Calibration: (V/MPa) <u>0.2767</u>		
62 mm V	LG: 19.73	HG: NA	Vert:		
<u>44 mm V</u>	LG: 32.3	HG: 309	Calibration: (V/MPa)		
22 mm V	Gain: 732.1				

Vessel Pressure:

Pore Fluid: _____

Calibrations (V/MPa)		Pressures (MPa)	Initial Voltage	Voltage @ Load
LG: 0.147	HG: 1.52	PpA:		
LG: 0.146	HG: 1.48	PpB:		
Gain : 0.1456		Pc:		
LG: NA	HG: NA	Pdiff:		

Data Logger Used: 8-channel

Control File

Horz. DCDT: ___ Long rod ___ Short rod

Vert. DCDT: ___ TT 2" Gain: High/Low
(HG: 0.57 mm/V LG: 2.85 mm/V)

(LR - HG: 0.622 mm/V LG: 1.27 mm/V

SR - HG: 0.64 mm/V LG: 1.32 mm/V)

Purpose/Description: _____

Acoustics blocks used L-Block

Temperature (°C): _____ Relative Humidity (%): _____

@ Hyd. Power Supply (HPS)

Chilled water at HPS

Chiller Unit

Process water at Chiller

14. Tank Temp (°C): _____

1. Temp In (°F): _____

6. Panel Temp (°F): _____

10. Temp In (°F): _____

15. Temp. Out (°C): _____

2. Pres. In (psi): _____

7. Panel Pres. (psi): _____

11. Pres. In (psi): _____

16. Pres. Out (psi): _____

3. Temp Out (°F): _____

8. Near Pres. In (psi): _____

12. Temp Out (°F): _____

4. Pres. Out (psi): _____

9. Near Pres. Out (psi): _____

13. Pres. Out (psi): _____

5. Flow (lpm): _____

3333 @ 5 MPa

3430 \uparrow 100 Hz, check Piezo-stack \rightarrow Good!

7600 \uparrow 10 kHz, Run 1

24197615 \downarrow 1 Hz, \uparrow NS, \uparrow 10 kHz, Run 2

30347700 \downarrow 1 Hz, \uparrow NS, \uparrow 10 kHz, Run 3

A2007790 \downarrow 1 Hz, \uparrow NS, \uparrow 100 Hz, Run 4 offset = 1.422

42377400 \downarrow 1 Hz, \uparrow 100 Hz, Run 5

42478900 \downarrow 1 Hz, \uparrow 100 Hz, Run 6

42609400 \downarrow 1 Hz, \downarrow NS \approx 1 MPa, Disp Mode, leave overnight

Piezo Stack Calibration - 3.75V
MPa

R	T
11	
8	8
5	5
2	2

PART B

1100 \uparrow NS = 10 MPa

2110 \uparrow 10 kHz, Run 7

2005160 \downarrow 1 Hz, \uparrow 10 kHz, Run 8

30052370 \downarrow 1 Hz, \uparrow 10 kHz, Run 9

42012400 \downarrow 1 Hz, \uparrow 100 Hz, Run 10

4237400 \downarrow 1 Hz, \uparrow 100 Hz, Run 11

4248600 \downarrow 1 Hz, \uparrow 100 Hz, Run 12

42612100 \downarrow 1 Hz

RST unplugged during last 10s of run? \rightarrow FIXED

PART C

108 @ 15 MPa

200 \uparrow 10 kHz, Run 13

20490200 \downarrow 1 Hz, \uparrow 10 kHz, Run 14

30800400 \downarrow 1 Hz, \uparrow 10 kHz, Run 15

42280600 \downarrow 1 Hz, Panel \rightarrow Ext 1, \uparrow 100 Hz, Run 16

42637802 \downarrow 1 Hz \uparrow 100 Hz, Run 17

42744751 \downarrow 1 Hz \uparrow 100 Hz, Run 18

42876625 \downarrow 1 Hz

42876880 \downarrow NS \approx 1 MPa, Disp Mode, lock, leave overnight

PART D

#520 @ 20 MPa

#600 10 kHz, RUN 19

#20000700 ↓ 1 Hz, 10 kHz, RUN 20

#30040800 ↓ 1 Hz, 10 kHz, RUN 21

#42091030 ↓ 1 Hz, 10 kHz, RUN 22

#42456200 ↓ 1 Hz, 100 Hz, RUN 23

#42561400 ↓ 1 Hz, 100 Hz, RUN 24

#42691700 ↓ 1 Hz, lock, Ext 1 → Panel, ↓ NS ≈ 1 MPa

Disp Mode, unplug Piezo-stack.

#42691920 Unload, End of Experiment

Oscillations Protocol

<i>Piezo-Stack Oscillations</i>
Amp1 = [0.2, 0.4, 0.6, 0.4, 0.8, 0.4, 1.0]MPa @ 10 Hz
Amp2 = [0.2, 0.4, 0.6, 0.4, 0.8, 0.4, 1.0]MPa @ 100 Hz
Amp3 = [0.2, 0.4, 0.6, 0.4, 0.8, 0.4, 1.0]MPa @ 200 Hz
F1 = [10, 50, 100, 200, 250]Hz @ 0.4 MPa
F2 = [10, 50, 100, 200, 250]Hz @ 1.0 MPa
<i>Biax Oscillations</i>
Amp4 = [0.2, 0.4, 0.6, 0.4, 0.8, 0.4, 1.0]MPa @ 0.1 Hz
Amp5 = [0.2, 0.4, 0.6, 0.4, 0.8, 0.4, 1.0]MPa @ 1.0 Hz
Amp6 = [0.2, 0.4, 0.6, 0.4, 0.8, 0.4, 1.0]MPa @ 10 Hz
F3 = [0.1, 1.0, 10]Hz @ 0.4 MPa
F4 = [0.1, 1.0, 10]Hz @ 1.0 MPa

<i>Experiment</i>			
Run1: ✓ Amp1, Amp2, Amp3	Run7: ✓ Amp1, Amp2, Amp3	Run13: ✓ Amp1, Amp2, Amp3	Run19: ✓ Amp1, Amp2, Amp3
Run2: ✓ F1, F2	Run8: ✓ F1, F2	Run14: ✓ F1, F2	Run20: ✓ F1, F2
Run3: ✓ Amp2, F2	Run9: ✓ Amp2, F2	Run15: ✓ Amp2, F2	Run21: ✓ Amp2, F2
Run4: ✓ Amp4, Amp5, Amp6	Run10: ✓ Amp4, Amp5, Amp6	Run16: ✓ Amp4, Amp5, Amp6	Run22: ✓ Amp4, Amp5, Amp6
Run5: ✓ F3, F4	Run11: ✓ F3, F4	Run17: ✓ F3, F4	Run23: ✓ F3, F4
Run6: ✓ Amp5, F4	Run12: ✓ Amp5, F4	Run18: ✓ Amp5, F4	Run24: ✓ Amp5, F4



Normal Stresses	
MPa	Volt
5	1.422
10	2.805
15	4.188
20	5.572


✓
✓
✓
✓

piezo stick calibration - 3.75 V/MPa

Internal cables


Surface

Vessel



Surface

Vessel



Acoustic cables

	60	
	59	
	58	
	57	

	8	
	7	
	6	